

Nutrition Information

		100 g	0,63 g (1 tablet)
• Energy			
	kJ	1352	9
	kcal	321	2
• Fat	g	3	<0,1
of which saturates	g	3	<0,1
• Carbohydrate	g	0	0
of which sugars	g	0	0
• Fibre	g	15	0,1
• Protein eqv.	g	67	0,4
of which			
L-Arginine	g	80	0,5
	mmol		2,9
• Salt	g	0	0



Arginin⁵⁰⁰ minis

- ✓ 230 tablets = 145 g
- ✓ 500 mg (2,9 mmol) L-Arginine/1 tablet



☎ 008000-9963829
(freecall from A, D, NL)
or + 49 (0) 8432 - 94 86 0

✉ service@metax.org
🛒 metax-shop.org
visit us on 🌐 metax.org

metaX Institut fuer Diätetik GmbH
Am Strassbach 5 · 61169 Friedberg/Germany



PRODUCT INFORMATION



FOOD FOR SPECIAL MEDICAL PURPOSES

Arginin⁵⁰⁰ minis

- for use in the dietary management of Urea Cycle Disorders or other medical conditions, when supplementation with L-Arginine is indicated

✓practical ✓flexible ✓discreet

- L-Arginine in the form of tablets
- neutral in taste and smell



✓ in case of e. g. Urea Cycle Disorders

- **Arginin⁵⁰⁰ minis** are an amino acid supplement in the form of coated tablets – exclusively with L-Arginine.
- **Arginin⁵⁰⁰ minis** contain 500 mg (2,9 mmol) L-Arginine per tablet.
- **Arginin⁵⁰⁰ minis** are a food for special medical purposes.
- **Arginin⁵⁰⁰ minis** are suitable for use in the dietary management of Urea Cycle Disorders, such as NAGS Deficiency, CPS Deficiency, OTC Deficiency, Citrullinemia Type I, adult-onset of Citrullinemia Type II or Argininosuccinic Aciduria, or other medical conditions, when supplementation with L-Arginine is indicated. Suitable from 3 years of age onwards.
- **Arginin⁵⁰⁰ minis** are not suitable in Arginase Deficiency (Argininemia).
- **Arginin⁵⁰⁰ minis** are highly concentrated.
- **Arginin⁵⁰⁰ minis** are not fortified with micronutrients.

✓ for use in the dietary management

- L-Arginine (**Arginin⁵⁰⁰ minis**) is used in the dietary management of Urea Cycle Disorders – when supplementation with L-Arginine is indicated.
- L-Arginine (**Arginin⁵⁰⁰ minis**) can also be used as a supplement in the dietary management of Propionic Aciduria with tendency to Hyperammonemia.
- In the dietary management of Pyridoxine-Dependent Epilepsy (PDE) with a Lysine restricted diet a supplementation with L-Arginine (**Arginin⁵⁰⁰ minis**) can also be indicated.
- **Arginin⁵⁰⁰ minis** can also be used in case of wound healing disorder.

✓ convenient and discreet

- **Arginin⁵⁰⁰ minis** are ready for application, no preparation necessary.
- **Arginin⁵⁰⁰ minis** can be consumed discreetly.

✓ neutral and flexible

- **Arginin⁵⁰⁰ minis** are neutral in taste and smell – they allow a tasteless intake of L-Arginine.

✓ easy and practical

- The surface of **Arginin⁵⁰⁰ minis** is very smooth. This is why you can easily swallow the tablets.
- For taking **Arginin⁵⁰⁰ minis** you simply need a glass of water or any other permitted drink: **Arginin⁵⁰⁰ minis** are practical, handy, easy!



✓ Dosage/Administration

- **Arginin⁵⁰⁰ minis** should always be taken with a sufficient quantity of liquid. They should be swallowed unchewed as whole tablets. They are not meant to dissolve in the mouth. We recommend to take the tablets in an upright position.
Note: The tablets are covered with a tasteless coating. This makes it easy to swallow them. However the layer is very thin. Thus, if the tablets were kept in the mouth for too long, this layer would dissolve and release the typical taste of Arginine. This is why we recommend the administration above.
- The daily dosage depends on age, body weight and individual medical condition and is determined under medical supervision.
- Also the right time for the daily intake should be determined by healthcare professionals.
- **Arginin⁵⁰⁰ minis** should at best be taken along with pre-calculated amounts of food.
- **Arginin⁵⁰⁰ minis** should be taken supplementary to the regular amino acid mixture – if applicable.
- The box with **Arginin⁵⁰⁰ minis** must be closed tightly after use, also with the inner lid.

✓ **Arginin⁵⁰⁰ minis**

✓ **Arginin⁵⁰⁰ minis**

✓ **Arginin⁵⁰⁰ minis**