



ZeroTP Piccos



In Short

- ZeroTP Piccos are a Food for Special Medical Purposes (FSMP) for use in the dietary management of Tyrosinemia or Alkaptonuria
- prolonged release of amino acids – Slow Release
- in form of very small tablets/pellets – coated
- suitable from 3 years of age
- highly concentrated protein supplement free of Tyrosine, with a small amount of Phenylalanine
- highly purified L-amino acids
- with vitamins, minerals and trace elements
- low in calories: no carbohydrate – small amount of fat
- very easy to swallow
- ready to use – ideal also for “ready to go”
- neutral in smell and taste
- in resealable and environmentally friendly plastic container (Securibox®), with dosing aid for individual and flexible dosing
- 1 dosing aid filled with 11,25 ml \cong 8,4 g product \cong 5 g protein eqv.
- 1 Securibox® contains 100 g ZeroTP Piccos \cong 60 g protein eqv.

Amino acids

ZeroTP Piccos	100 g	8,4 g
L-Alanine	9	5,1
L-Arginine	9	3,0
L-Aspartic acid	9	8,7
L-Cystine	9	1,6
L-Glutamine	9	2,3
L-Glutamic acid	9	6,3
Glycine	9	5,3
L-Histidine	9	1,6
L-Isoleucine	9	3,7
L-Leucine	9	5,9
L-Lysine	9	4,9
L-Methionine	9	1,7
L-Phenylalanine	mg	100
L-Proline	9	6,7
L-Serine	9	3,6
L-Threonine	9	3,7
L-Tryptophan	9	1,5
L-Tyrosine	9	–
L-Valine	9	5,5

Product Profile

ZeroTP Piccos are a highly concentrated mixture of highly purified L-amino acids free of Tyrosine, with a small amount of Phenylalanine, in form of very small tablets/pellets.

The amino acid mixture is embedded in sodium alginate. This leads to a slow/prolonged release of amino acids within the gastrointestinal tract – Slow Release.

Good to know Slow Release

The individual amino acids are gradually released within the gastrointestinal tract over an extended period of time. This contributes to a digestion of the amino acids similar to that of natural protein.

(see p. 2)

ZeroTP Piccos do not contain carbohydrate and only a small amount of fat, thus low in calories.

ZeroTP Piccos are a pressed and coated protein supplement.

ZeroTP Piccos allow a tasteless intake of the amino acid mixture. They do not leave an unpleasant aftertaste and do not cause breath malodour or belching.

The very small tablets/pellets can be consumed discreetly and directly, no preparation necessary.

ZeroTP Piccos offer more flexibility, e.g.

concerning intake outside the home and with regard to a wider choice of drinks suitable for intake, without drawing attention to the diet – a big plus with regard to good compliance.

ZeroTP Piccos provide an alternative for all those patients who wish to benefit from the advantages of our ZeroTP minis but cannot swallow them.

Since packed in a plastic container (Securibox®) with a dosing aid for individual and flexible measuring, handling is particularly easy: 1 dosing aid filled with 11,25 ml

\cong 8,4 g product \cong 5 g protein equivalent.

1 Securibox® contains 100 g ZeroTP Piccos \cong 60 g protein equivalent.

ZeroTP Piccos are supplemented with vitamins, minerals and trace elements.

Suitable from 3 years of age.

Phenylalanine Enrichment

ZeroTP Piccos contain a small quantity of phenylalanine, an essential amino acid required for protein biosynthesis, and thus for growth and hormone production (e.g., dopamine and adrenaline). Insufficient levels may result in neurological symptoms. Phenylalanine supplementation is recommended for patients with tyrosinaemia in order to:

- ensure an adequate supply of this essential amino acid,

- prevent deficiencies resulting from the low-phenylalanine diet needed,
- optimize metabolic balance.

The phenylalanine content has been adjusted to align with GMP-based products, an approach that has proven effective in routine clinical practice.

Convenience Product

ZeroTP Piccos are a convenient, ready to use option for the intake of an amino acid mixture, also when on the way: discreetly, in tablet form/pellets.

Administration


ZeroTP Piccos should always be taken with a sufficient quantity of liquid.

Measuring or Weighing?

The following applies to the enclosed dosing aid:

5 g protein eqv. \cong 11,25 ml

This corresponds to a filling between the 10 ml indication and the mark above it.



It might be helpful to mark the corresponding spot with a permanent marker.

Nevertheless, we recommend weighing the desired amount of Piccos if possible:

5 g protein eqv. \cong 8,4 g product.

Tip: If you check the first measurements using a scale (display with at least 1 digit after the decimal point), you will quickly get a feeling for how full the dosing aid is when the target value of 8,4 g is reached.

Intake

The tablets/pellets should be swallowed unchewed as whole. **ZeroTP Piccos** should always be taken with a sufficient quantity of liquid. The tablets/pellets are not intended to be dissolved in water or other drinks prior to intake. **ZeroTP Piccos** are not meant to dissolve in the mouth.

We recommend to take the tablets/pellets in an upright position.

ZeroTP Piccos should be used within three weeks – the original Securibox® being closed tightly after every use.

ZeroTP Piccos are coated with a flavour- and odour-neutral layer – therefore they do not produce an unpleasant aftertaste or breath malodour and there is no annoying belching.

Thanks to their protective coating the tablets/pellets are easy to swallow. This coating is extremely thin and water-soluble. If it were to dissolve, the tablets/pellets would taste like a typical amino acid mixture.

ZeroTP Piccos are tooth-friendly. Intake of the coated tablets/pellets with water (including mineral water) contributes to oral hygiene and helps prevent possible tooth damage caused by caries.

Environmentally friendly Packaging

ZeroTP Piccos are packed in a plastic container. The plastic of this so called Securibox® consists of renewable, plant-based raw material obtained from controlled, eco-friendly cultivation – sustainable and CO₂ efficient.

Function **ZeroTP Piccos** substitute that part of the protein in the diet which may not be taken up from natural food sources.

Indication **ZeroTP Piccos** are a food for special medical purposes and as protein supplement suitable for the dietary management of Tyrosinemia or Alkaptonuria.

Dosage The daily total amount of amino acid mixture depends on age, body weight and individual medical condition and should be re-examined and adjusted in accordance with the results of regular monitoring. The daily dosage of amino acid mixture should at best be divided into 3 – 5 single portions.

ZeroTP Piccos can be combined with other products from the ZeroTP-system. The diet must be supplemented with energy, natural protein and other nutrients in prescribed quantities.

Important Notice Must only be used under medical supervision. Not for use as a sole source of nutrition. For enteral use only. Only for people with proven Tyrosinemia or Alkaptonuria. Suitable from 3 years of age.

Reference:

- de Jesus Raposo MF, de Moraes AM, de Moraes RM. Emergent Sources of Prebiotics: Seaweeds and Microalgae. *Mar Drugs*. 2016.

GOOD TO KNOW Slow Release

Slow Release protein substitutes are characterised by a prolonged release of amino acids.

By embedding the amino acid formulation in alginate, a natural polymer from brown algae forming a hydrophilic gel, the amino acids are gradually released within the gastrointestinal tract over an extended period of time. This contributes to a digestion of the amino acids similar to that of natural protein. The prolonged absorption kinetics of Slow Release protein substitutes promote endogenous protein synthesis and prevent catabolism. In animal and human studies, Slow Release amino acid formulations have been shown to lead to improved metabolic control and lower peak concentrations of individual amino acids in the blood. Alginates exhibit a prebiotic effect by promoting the microbiome. In addition, various clinical studies have demonstrated that alginates have the potential to reduce body weight and cholesterol levels and to prevent diabetes by reducing post-prandial glucose and insulin peaks.

NUTRITION INFORMATION

ZeroTP Piccos	100 g	8,4 g
Energy		
kJ	1165	98
kcal	276	23
Fat		
g	1	0,08
of which saturates	g	1
Carbohydrate	g	0
of which sugars	g	0
Fibre		
g	13	1,1
Protein eqv.		
g	60	5
Amino acids		
g	72	6
Salt		
g	2,4	0,2
Vitamins		
		g protein
Vitamin A	µg	1080
		91
		18
Vitamin D3	µg	39
		3
		0,6
Vitamin E	mg	21
		1,8
		0,35
Vitamin K1	µg	77
		7
		1,3
Vitamin C	mg	126
		11
		2
Thiamin (Vit. B1)	mg	2
		0,17
		0,035
Riboflavin (Vit. B2)	mg	2,4
		0,2
		0,04
Niacin	mg	48
		4
		0,8
Vitamin B6	mg	2
		0,16
		0,03
Folic acid	µg	360
		30
		6
Vitamin B12	µg	6
		0,5
		0,1
Biotin	µg	48
		4
		0,8
Pantothenic acid	mg	9
		0,8
		0,15
Minerals		
Sodium	mg	950
		80
		16
Potassium	mg	1514
		127
		25
Calcium	mg	1507
		127
		25
Phosphorus	mg	976
		82
		16
Magnesium	mg	372
		31
		6
Trace elements		
Iron	mg	22
		1,8
		0,36
Zinc	mg	22
		1,8
		0,36
Copper	mg	2,4
		0,2
		0,04
Manganese	mg	7
		0,6
		0,11
Fluroide	mg	0,9
		0,08
		0,015
Selenium	µg	87
		7
		1,5
Chromium	µg	120
		10
		2
Molybdenium	µg	102
		9
		1,7
Iodine	µg	270
		23
		4,5
FURTHER NUTRITION INFORMATION		
L-Carnitine	mg	54
		5
		0,9
Choline	mg	540
		45
		9
myo-Inositol	mg	282
		24
		5
Taurine	mg	60
		5
		1

INGREDIENTS

Bulking agents: sodium alginate & microcrystalline cellulose & silicium dioxide & isomalt, L-lysine-L-aspartate, potassium-L-glutamate, L-proline, L-leucine, L-valine, glycine, L-arginine-L-aspartate, L-alanine, calcium phosphate, L-isoleucine, L-threonine, L-serine, magnesium L-aspartate, L-glutamine, N-acetyl-L-methionine, anti-caking agents: calcium silicate & magnesium stearate, glazing agents: E 464 & E 463 & E 422, colouring food: safflower & lemon, colour: copper complexes of chlorophylls, L-cystine, L-histidine, L-tryptophan, calcium carbonate, L-phenylalanine, choline, L-carnitine, vitamin E, vitamin A, sodium molybdate, folic acid, inositol, niacin, ferrous sulphate, sodium selenite, vitamin C, zinc sulphate, pantothenic acid, taurine, vitamin K, sodium fluoride, manganese sulphate, chromium (III) chloride, vitamin B2, vitamin B1, vitamin B6, potassium iodide, cupric sulphate, biotin, vitamin D, vitamin B12.

Delivery Unit	box, Securibox® 7 x Securibox® at 100 g = 700 g
Article Number	xx-001-31070
Delivery to	Pharmacies, clinics
Storage	Store in a cool, dry place.



metaX Institut für Diätetik GmbH
Am Strassbach 5
61169 Friedberg/Germany

phone +49 (0) 8432 - 94 86 0
e-mail service@metax.org
web metax-shop.org · metax.org